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Basic Research  
EA02010  
poster presentation

## Correlation of drilling frequency with torque, and with vertical force ドリル回転数とトルク、垂直的力の関係

Yoshihiro Sakaguchi, Masashi Takahashi GC Corporation, Tokyo, Japan

### Background and Aim

There are various factors for success of the implant treatment. One of the keys is to make an appropriate hole by using an appropriate drill. Each implant system has its recommended drill frequency which is unique to each system. However it is not clear how the drill frequency should be set.

The aim of this study is to evaluate the correlation of drilling frequency with torque, and with vertical force.

### Methods and Materials

The maximum torque (Max Mz), maximum vertical force (Max Fz) and minimum vertical force (Min Fz) was measured by 4-Component Dynamometer (Kistler Japan Co., Ltd. Fig.1) by drilling up Polytetrafluoroethylene block by the diameter 2.0mm twist drill.

Frequencies of drill were 500rpm, 1000rpm, 1500rpm and 2000rpm. The measurement was performed three times by each rotational speed.

4-Component Dynamometer (Kistler Japan Co., Ltd. Fig.1) Multiple quartz plates can be connected electrically and mechanically in parallel. 4-component sensor is created for Fx, Fy, Fz and Mz.


Fig.1 4-Component Dynamometer

Table 1 Time base range and Threshold level of 4-Component Dynamometer

Time base range	Horizontal force (FX Fy)	-5~5KN
	Vertical force (Fz)	-5~20KN
	Torque (Mz)	-200~200N·m
Threshold level	Horizontal force (FX Fy)	<0.01N
	Vertical force (Fz)	<0.02N
	Torque (Mz)	<0.02N·cm

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
Fig. 2 Correlation of drilling frequency with torque


Fig. 3 Correlation of drilling frequency with vertical force

### Conclusions

The cutting torque increases when the frequency increases. Maximum and minimum vertical force decreases when the frequency increases. These results show If frequency is too high, we should be careful for the bone burn. If frequency is too high or too low, we should be careful for the over drilling.

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- 1) TG Wilson Jr. J Periodontol 2009: The positive relationship between excess cement and peri-implant disease: a prospective clinical endoscopic study

Basic Research  
A02014  
poster presentation

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Tsutom Abe, Masaya Seto, Takaharu Takeshita GC Corporation, Tokyo, Japan

## Introduction

The abutments for Dental Implant system are classified into two types. One is the screw-retained type (Fig.1) and the other is the cement-retained type (Fig.2). Especially, screw-retained abutments with superstructure (prosthetic screw and gold cylinder with casting alloy) have a complicated structure compared with the cement-retained abutments, so there are a few reports of the fatigue strength of those abutments. In this study, the fatigue strength of the implant system including implant, screw-retained abutment and superstructure with several angulations is investigated.

Fig.1 Screw-retained abutment  
Fig.2 Cement-retained abutment

## Material and Methods

### Materials

#### Fatigue Test Setting and Materials

①	Components	Materials
②	Implant	ASTM F136
③	Abutment	ASTM F136
④	Gold Cylinder	Gold Alloy
⑤	Prosthetic Screw	Gold Alloy
⑥	Casting Alloy	Gold Alloy (Type 4)
⑦	Cap	Cr-Mo Steel

#### Components of Trial Implant system

##### 1) Superstructure

A) Wax-up models were made on the Gold Cylinder  
B) Substituted by Casting Alloy and heat-treated of the porcelain baking process.

The shape of superstructure was the mechanically worst case determined by Finite Element Method (FEM), as below.

##### 2) Screw-retained Abutments

A) Straight type Abutment  
B) 30° Angled Abutment

##### 3) Implant

Two shapes of Abutments were chosen


Trial implants have an Internal butt joint connection

## Experimental procedure

### FEM


The shape of superstructure was the mechanically worst case determined by Finite Element Method (FEM).

- Software: Solidworks Simulation 2011 SP5.0
- Simulated Shapes: a) Straight shape (Long), b) Straight shape (Short), c) 10° Angled shape (Long)
- Loading: 200N


### Fatigue Test

- In accordance with ISO 14801
- Machine: MMT-500NV-10 (Shimadzu)
- Settings of Angulation:
  - 30° angled (straight abutment),
  - 50° angled (Angled abutment)
- Frequency: 15Hz
- Loading Cycles: 5 million cycles
- Maximum endured load was determined by three specimens reached the specified number of cycles with no failures.

## Results

### FEM

- The position of stress concentration: The root of Abutment thread was indicated in all shapes.
- Maximum stress: Similar stress was shown in all shapes at 200N loaded.
- Selected shapes of fatigue tests
  - Straight shaped (Short) superstructure for Straight Abutment
  - Angle shaped superstructure for Angled Abutment

### Fatigue Test



Fig.10 Locations, Fractured positions, and load-cycle diagrams

The fatigue tests of angle type were started from 300N which was determined by flexural test. Since there were no fractures at the 300N, the fatigue tests were performed at higher loads. As a result, maximum endured load was 375 N with no fractures. At the 400N, the fractures of implant and abutment screw were found. The fractured position was consistent with the position of stress concentration indicated by FEM. On the other hand, the fatigue tests of straight type were started from 375N maximum of angle type. At 375N, the fracture of abutment screw was occurred, so the fatigue tests were performed at lower loads. As a result, maximum endured load was 230N with no fractures.

## Conclusions

Generally, implants with greater insertion angle are more likely to be suffered from high flexural stress. According to the result of this study, the maximum endured load of angle type abutment was shown to be higher than that of straight type abutment. Therefore, it was indicated that load to implant systems can be controlled by adjusting direction with appropriate selection of the angle type abutment.

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