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●Reference

- Misch CE. Contemporary implant dentistry. 3ra. Edición. St. Louis, Missouri: Mosby Elsevier; 2008.
- Gonzales J. The evolution of dental materials for hybrid prosthesis. The Open Dentistry Journal. 2014; 8: 85-94.
- Zarb GA. The longitudinal clinical efficacy of osseointegrated implants a 3 year report, International Journal Oral and Maxillofacial Implants. 1987; 2: 91-100.
- Cobb G, Metcalf M, Parsell D. An alternative treatment method for a fixed-detachable prosthesis. A clinical report. The Journal of Prosthetic Dentistry. 2003; 89 (3): 239-243.
- Shibli JA, Piatelli A, Lezzi G. Effect of smoking on early bone healing around oxidized surfaces: a prospective, controlled study in humans jaw. J Periodontol. 2010; 81:575-583.
- Bain C, Moy P. The association between the failure of dental implants and smoking. Int J Oral Maxillofac Implants. 1993; 8: 609-615.
- Grunder U, Gaberthuel T, Boitel N. Evaluating the clinical performance of the osseotite implant: Defining prosthetic predictability. Compend Contin Educ Dent. 1999; 20: 628-640.
- Bain C. Long-term satisfaction in dental implant patients. Toronto: 16th Scientific Meeting Academy of Osseointegration; 2001.
- Balarezo JA. Prótesis sobreimplantes en el edéntulo total: Planificación y elaboración. Lima, Perú: Editorial Savia. 2014;.
- Brida A, Agar J. A classification system of patients for esthetic fixed implant-supported prostheses in the edentulous maxilla. Compend Contin Educ Dent. 2010;31 (5):366-8.
- Tarnow DP, Emtiaz S, Classi A. Immediate loading of threaded implants at stage 1 surgery in edentulous arches: Ten consecutive case reports with 1-to 5-year data. Int J Oral Maxillofac Implants. 1997;12(3):319-2
- Real-Osuna J, Almendros-Marqués N, Gay- Escoda C. Prevalence of complications after the oral rehabilitation with implant-supported hybrid prostheses. Med Oral Patol Oral Cir Bucal. 2012;17 (1): 16-21.
- Qamheya AH, Yenyol S, Arisan V. Full Mouth Oral Rehabilitation by Maxillary Implant Supported Hybrid Denture Employing a Fiber Reinforced Material Instead of Conventional PMMA. Case Rep Dent. 2015;2015:841745
- Kwon T, Bain P, Levin L. Systematic review of short (5-10 years) and long-term (10 years or more) survival and success of full-arch fixed dental hybrid prostheses and supporting implants. Journal of Dentistry, 2014; 42 (10): 1228-1241.
- Egilmez F, Ergun G, Cekic-Nagas I, Bozkaya S. Implant-supported hybrid prosthesis: Conventional treatment method for borderline cases. Eur J Dent. 2015; 9: 442-8.

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略歴・所属団体©1997～2002年 マドリード欧州大学。2004～2006年 マドリード欧州大学歯科インプラント学修了。2005～2007年 マドリード・コンプルテンセ大学口腔生物学修士。2005年 Multidisciplinary Dental Practice CIMA 設立 (スペイン、マドリード)。2008年 マドリード・コンプルテンセ大学審美歯科卒業プログラム講師。2009年 マドリード欧州大学歯周病学卒業プログラム准教授。2014年 FORプラットフォーム専門家