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Q1 Despite the ongoing evolution of bonding agents still many dentists prefer to use either one or the other type of the older generation of multi-step adhesives. What do you think is the reason for that?

Bart Many dentists still prefer a certain type of bonding procedure, and feel familiar and are satisfied with a specific product; therefore, they do not see a clear need to change their bonding procedure and dental adhesive. In my opinion, what they were taught in their studies also contributes to their choice. Dentists tend to continue doing it that way if they believe the technique they have learned is efficient.

Q2 For a long time the scope of research was focused on the variety of adhesive monomers and the methodology of bond strength testing as well as the level of bond strengths itself. Nowadays the scope seems to be more on bond durability and long-term performance. How do you consider the different generations of bonding agents on that point?

Bart Most adhesives nowadays guarantee sufficient 'immediate' bonding effectiveness, which means they perform well when bonded to model dentin under ideal laboratory circumstances and tested relatively shortly after their application, usually after 24 hours or 1 week. Consequently, in most cases we cannot differentiate among today's adhesives based only on immediate bond-strength data.

For this reason, it is extremely important to challenge the products to learn their behavior in the long term. As editor of the *Journal of Adhesive Dentistry*, I must tell you that we do not accept papers reporting bond-strength data without aging. Of course, artificial aging in the lab cannot be directly translated to the clinical situation but, at least, it gives an idea on the stability of the bond.

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
Though a scientifically documented product-dependency exists, the 3-step ER and 2-step SE adhesives are still considered the gold standards. This is mainly because sufficient and repeatedly consistent evidence of favorable long-term bonding performance was recorded for particular commercial adhesives in both laboratory and clinical research, the latter beyond ten years of clinical service. All other simplified adhesives should somewhat be regarded as 'trade-off' adhesives, which traded in some of their bonding performance for ease-of-use, hereby not yet reaching the proven level of the gold standards, while quite often also not having been tested clinically for a sufficiently long time.

Q3 There are perceptions about a potential weakness of simplified E&R adhesives in regard to collagen degradation and for both strategies, E&R and SE approaches, in regard to the hydrolysis of the polymer matrix itself. What is important for the general practitioner to know on these topics?

Bart Two main bond-degradation pathways have been documented in scientific literature. There is first the so-called enzymatic bio-degradation (1), involving MMPs along with cysteine cathepsins. A lot of attention in research is given to this pathway; however, the proportional contribution of such enzymatic activity to bond degradation is still unclear and in my personal opinion still much less relevant than the second bond-degradation pathway. Second, there is water sorption (2), or water uptake by the adhesive interface, inducing hydrolytic bond-degradation mechanisms.

I definitely consider water sorption as the most important degradation pathway. The oral environment and the dentin itself are the two water sources. The adhesive interface can absorb water from the underlying dentinal tissue through osmosis, and this may lead to hydrolytic degradation of the adhesive interface with dentin.

I believe that nowadays the importance of a good interfacial seal to prevent water sorption is underestimated, while this is besides the actual bonding performance another primary requirement an adhesive should meet. Too much attention goes nowadays to enzymatic bio-degradation, which I consider a kind of research 'trend' that will gradually pass with time.

Most important for bond durability is a stable and hydrophobic adhesive interface that ensures a good seal and prevents/limits water uptake. It must be clear for the general practitioners that adhesives should provide this hardly permeable seal.

Q4 Tooth erosions, bulk-fill techniques, shorter light-curing protocols, just to mention some trends, implicate new challenges for adhesive dentistry. What are your recommendations in order to achieve a high level of bond strength, durability and long-term performance in direct protocols?

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Bart My preferred bonding protocol combines selective enamel etching with a SE approach based on 10-MDP. Among many functional monomers, 10-MDP is nowadays considered the most effective monomer. On top of micromechanical interlocking achieved by 10-MDP's etching capacity, this monomer provides additional chemical interaction potential. Chemical interaction is the innermost contact that molecules can have to each other and that's fundamental for bond durability. A common misunderstanding is the idea that chemical interaction will increase bond strength; it will not! Chemical interaction is responsible for maintaining bond strength at a stable level.

Also of high importance for bond-strength durability is the immediate seal of the interface. Postponing light curing creates osmotic effects with the consequence that water is absorbed from the underlying dentin. As mentioned earlier, this targeted interfacial seal must be, of course, highly hydrophobic.

I also strongly believe that in order to protect and stabilize the hybrid layer, the adhesive layer should have a certain thickness. A thick adhesive layer can act as a shock-absorbing layer or stress absorber, which for instance compensates for the shrinkage of the overlaid composite and so will impose tensile stress to the underlying adhesive interface.

In summary, I would like to state that selective enamel etching, a SE approach based on 10-MDP, a hydrophobic seal and an adhesive layer with a certain thickness are essential requirements to achieve long-term bond stability.

Q5 GC has a successful history in the development of HEMA-free bonding agents. Where do you see the specific benefits of this feature?

Bart Today, HEMA is still added to many commercial adhesives due to some positive features. It can act as a co-solvent for other monomers in preventing water/monomer phase separation. Thanks to its low molecular weight and thus small size, along with its high hydrophilicity, HEMA is an effective surface-wetting as well as interdiffusion/penetration agent, the latter to infiltrate into the wet, demineralized collagen-rich dentin relatively deeply exposed upon phosphoric-acid etching as part of an E&R bonding protocol.

However, HEMA has also major disadvantages. It is too hydrophilic and this leads to high water sorption. Many commercial adhesives contain water and HEMA in the same solution, and the fact that HEMA is too compatible with water makes it nearly impossible to remove water from the interface. Another negative point is HEMA's low polymerization ability, remaining soft as a kind of gel despite being cured. If it does not cure properly, the resin matrix will consequently not be strong; so, HEMA has a low contribution to the mechanical strength and stability of the polymer matrix, and actually weakens it. HEMA is also very sensitive to hydrolysis. Another negative aspect is that, as it is a small molecule, it can compete with other functional monomers, inhibiting chemical interaction; for example, our research has shown that HEMA can affect 10-MDP's nano-layering. Its unfavorable biocompatibility, particularly in terms of its documented allergic potential, is finally another disadvantage of HEMA that cannot be ignored.

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
※10- methacryloyloxydecyl dihydrogen phosphate

Due to all these limitations, dental manufacturers attempt to develop and produce HEMA-free/poor adhesives, by which there is definitely a need to find an effective replacement for HEMA.

Q6 What is your wishlist in regard to a most effective 2-step bonding agent in terms of composition, technological features and performance?

Bart A separate primer should serve as adhesion promoter. This primer has to be 10-MDP based, rendering chemical binding potential. This 10-MDP monomer must be of good quality and be present in the solution in a sufficiently high concentration. The primer has to contain also photo-initiators, not enabling the primer to be cured of course, but to make sure that all areas, even in the deeper parts of the hybrid layer, will receive photo-initiators, so the infiltrated resin can polymerize efficiently in situ after the adhesive resin in the third and final step is applied.

The final step involves indeed the copious application of a separate adhesive resin, the actual glue. This adhesive should provide a good seal of the interface, it has to polymerize well, must present good mechanical properties and must be of hydrophobic nature to reduce water uptake. And it must be applied in a sufficiently thick layer.

We believe, based on previous researches, that a thicker adhesive layer can serve as a stress/shock-absorbing layer. When the composite that is placed on top of this layer shrinks, the interface is subjected to tensile stress; with an elastic layer (buffer layer) we can prevent that this tensile stress will cause interface de-bonding.

Q7 In your opinion, the application of an HEMA-free hydrophobic layer following the application of an hydrophilic primer would optimize the adhesive performance to the tooth substrate ?

Bart Due to the water content, dentin must always be properly sealed to avoid that water is absorbed by the adhesive interface through osmosis, which causes interfacial porosities and consequently weakens the interface.

This sealing has to be realized immediately and the application of a HEMA-free hydrophobic layer, when correctly polymerized, can reduce water uptake and hydrolytic bond degradation.

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
Message from Professor

G2-Bond Universal (GC), a universal adhesive in a new configuration for stable long-term bonding to tooth enamel and dentin

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(BIOMAT – Biomaterials Research group, Department of Oral Health Science, KU Leuven)

To durably bond to tooth tissue, believed primary important is to achieve an interface of highly hydrophobic nature that keeps the adhesive interface sealed against water uptake. Water sorption and resultant hydrolytic effects constitute the main bond-degradation pathway. In this regard, the current generation of universal adhesives (UAs) should be regarded as 'trade-off' adhesives; they combine the primer and adhesive-resin function in a single solution that proportionally contains more solvent and less resin, as compared to when a separate primer and adhesive resin are employed by multi-step adhesives (Van Meerbeek et al., 2020). UAs typically present with a thin film thickness without stress-absorbing potential (Ahmed et al., 2020). Having traded in some bonding performance for ease-of-use, UAs in their current configuration have not yet approached the gold-standard multi-step adhesive benchmarks. A recent multi-parameter research project revealed that the experimental precursor of G2-Bond Universal, instructing to separately apply a 10-MDP-based primer prior to a highly hydrophobic particle-filled adhesive resin in both etch-and-rinse or self-etch mode, successfully combined high bonding potential with bond-degradation resistance thanks to limited water uptake of a significantly thick adhesive film (Yao et al., 2020). This new UA configuration promises durable clinical adhesive performance.

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Van Meerbeek B, Yoshihara K, Van Landuyt K, Yoshida Y, Peumans M. From Buonocore's Pioneering Acid-Etch Technique to Self-Adhering Restoratives. A Status Perspective of Rapidly Advancing Dental Adhesive Technology. *J Adhes Dent.* 2020;22(1):7-34.

Yao C, Ahmed MH, Li X, Nedeljkovic I, Vandooren J, Mercelis B, Zhang F, Van Landuyt KL, Huang C, Van Meerbeek B. Zinc-Calcium-Fluoride Bioglass-Based Innovative Multifunctional Dental Adhesive with Thick Adhesive Resin Film Thickness. *ACS Appl Mater Interfaces.* 2020;12(27):30120-30135.